




CLOSE

Refractometer MASTER-M



Shop



Send an inquiry



Quotation Request

Can be widely used for beverages, fruit juices, soups, beef bowl sauce, water-soluble cutting oil, etc. MASTER-M is made of metal, and MASTER-PM is plastic. For salty and acidic samples, choose the plastic body, MASTER-PM. This product is the cheapest type with no temperature compensation function. Since the concentration changes as the measurement temperature changes, MASTER- α and MASTER-T are recommended if you have trouble managing the temperature at each measurement.

Specifications

Model	MASTER-M
Cat.No.	2313
Scale	Brix
Range	Brix : 0.0 to 33.0%
Minimum Scale	Brix : 0.2%
Accuracy	Brix : $\pm 0.2\%$
Repeatability	Brix : $\pm 0.1\%$
Dimensions & Weight	3.2 x 3.4 x 20.3cm, 155g

Options

- SUPER EXTRACTOR : RE-29401
- 10% Sucrose ($\pm 0.03\%$) : RE-110010
- 20% Sucrose ($\pm 0.03\%$) : RE-110020
- 30% Sucrose ($\pm 0.03\%$) : RE-110030
- Daylight plate for small volume samplings : RE-2311-67M
- Daylight plate : RE-2315-60M
- * Calibration Certificate : Contact an ATAGO representative for details.

